

Horus HD Digital Scope

Multi integrated imaging system.

The Horus HD Digital Scope System is an easy-to-use hand-held video system used for capturing images of the body. This system utilises interchangeable attachments making it ideal for multi discipline medical applications.

Attachments for ear nose & throat, dermatology ophthalmology, audiology, women's health and general examinations are just a few of the potentials of this unique system.

ear nose & throat
ophthalmology
general viewing
dermatology



Horus HD Digital Scope



1080p High Definition images are captured and displayed on the 3.5" / 90mm full colour LCD display.



Illumination is provided by powerful integrated LEDs.



The fine focus control is conveniently located where your index finger comfortably rests on the handle.



Conveniently located controls make this system extremely user friendly. Most functions can be controlled with just your thumb.



Each specialty attachment can be quickly and easily changed by just twisting the optic and removing it from the scope handle. Go from otoscopy to dermatology in just seconds.

Specifications: Series 1

Resolution	1920 x 1080p HD
Light Source	Light emitting diode (LED)
Display	3.5" / 90mm full colour LCD
Image Format	Jpeg (still) & H.264 (video)
Media Card	Micro SD 2GB (included)
Power Source	Rechargeable lithium-ion battery
Operating Time	3 hours
Weight	340g / 12oz
Focusing Range	5mm to 30mm

Specifications: Series 3

Resolution	2560 x 1920p HD
Light Source	Light emitting diode (LED)
Display	3.5" / TFT LCD
Image Format	Jpeg (still) & H.264 (video)
Media Card	SD 8GB (included)
Power Source	Rechargeable lithium-ion battery
Weight	275 grams
Focusing Range	7-20mm

For further information and a demonstration, contact **Visionflex** at info@visionflex.com.au

1300 059 926 +61 2 8438 9898

visionflex.com.au

visionflex
Equality healthcare